



LARYNGOSCOPE
ENT DEVICES
CATALOGUE



PS-L-001

Handle Standard Illumination
 Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ 20-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 20-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 20-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 20-200-41

Sets available in hard case or soft pouch



20-100-00

Mc - Intosh Blade Standard Illumination
 Size: 0



Mc - Intosh Blade Standard Illumination
 Size: 1



Mc - Intosh Blade Standard Illumination
 Size: 2



Mc - Intosh Blade Standard Illumination
 Size: 3



Mc - Intosh Blade Standard Illumination
 Size: 4

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Mc – Intosh Laryngoscopes
Fiber Optic Illumination
 Adult Size



Mc - Intosh Blade Fiber Optic Illumination
 Size: 0



Mc - Intosh Blade Fiber Optic Illumination
 Size: 1



Mc - Intosh Blade Fiber Optic Illumination
 Size: 2



Mc - Intosh Blade Fiber Optic Illumination
 Size: 3



Mc - Intosh Blade Fiber Optic Illumination
 Size: 4

PS-L-002

Handle Fiber Optic Illumination
 Size: Adult
 Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ 21-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 21-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 21-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 21-200-41

Sets available in hard case or soft pouch

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



PS-L-003

Handle Standard Illumination
 Size: Child / Infant

Handle options: Adult & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Child Handle	▪ 22-200-30
▪ Set of 3 Blades Size 0, 1, 2 and Child Handle	▪ 22-200-31
▪ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	▪ 22-200-40
▪ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	▪ 22-200-41

Sets available in hard case or soft pouch



Miller Blade Standard Illumination
 Size: 00



Miller Blade Standard Illumination
 Size: 0



Miller Blade Standard Illumination
 Size: 1



Miller Blade Standard Illumination
 Size: 2



Miller Blade Standard Illumination
 Size: 3

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Miller Laryngoscopes
Fiber Optic Illumination
 Adult Size



Miller Blade Fiber Optic Illumination
 Size: 0



Miller Blade Fiber Optic Illumination
 Size: 1



Miller Blade Fiber Optic Illumination
 Size: 2



Miller Blade Fiber Optic Illumination
 Size: 3



Miller Blade Fiber Optic Illumination
 Size: 4

PS-L-004

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ 23-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 23-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 23-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 23-200-41

Sets available in hard case or soft pouch

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Mc – Intosh Integrated Laryngoscopes

Fiber Optic Illumination

Adult Size



PS-L-005

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Adult Handle	▪ 24-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 24-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 24-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 24-200-41

Sets available in hard case or soft pouch



Mc - Intosh Blade Fiber Optic Illumination
Size: 00



Mc - Intosh Blade Fiber Optic Illumination
Size: 0



Mc - Intosh Blade Fiber Optic Illumination
Size: 1



Mc - Intosh Blade Fiber Optic Illumination
Size: 2



Mc - Intosh Blade Fiber Optic Illumination
Size: 3



24-100-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

McCoy Flexi - Tip Laryngoscopes

Fiber Optic Illumination

Adult Size



PS-L-006

Handle Fiber Optic Illumination
Size: Adult

Standard Sets

Contents	Code
▪ Set of 2 Blades Size 2, 3 and Adult Handle	▪ 25-200-30
▪ Set of 2 Blades Size 3, 4 and Adult Handle	▪ 25-200-31
▪ Set of 3 Blades Size 2, 3, 4 and Adult Handle	▪ 25-200-40

Sets available in hard case or soft pouch



McCoy Blade Fiber Optic Illumination
Size: 1



McCoy Blade Fiber Optic Illumination
Size: 2



McCoy Blade Fiber Optic Illumination
Size: 3



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



PS-L-007
Handle Fiber Optic Illumination
Size: Adult



McCoy Blade Fiber Optic Illumination
Size: 1



McCoy Blade Fiber Optic Illumination
Size: 2



McCoy Blade Fiber Optic Illumination
Size: 3



Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	26-200-30
Set of 2 Blades Size 3, 4 and Adult Handle	26-200-31
Set of 3 Blades Size 2, 3, 4 and Adult Handle	26-200-40

Sets available in hard case or soft pouch

Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



PS-L-008

Brass Otoscope Fiber Optic
Brass Otoscope Conventional



PS-L-009

ENT Set Complete



PS-L-010

Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available
- Fiber Optic Light Source



PS-L-011

Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



Color Options

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source

Standard Accessories

Nasal Speculum



PS-L-012



PS-L-013



PS-L-014

Bulbs



PS-L-015



PS-L-016



PS-L-017



PS-L-018

Ear Speculas



PS-L-019



PS-L-020



PS-L-021



PS-L-022

Ear Speculas & Testing Tube



PS-L-023



PS-L-024



PS-L-025



PS-L-026



PS-L-027



PS-L-028



PS-L-029