

PENTAX SURGICAL

Manufacturer of all kind of Laryngoscope and ENT Devices Instruments

2026-2030 CATALOGUE | EDITION



A photograph of a modern, multi-story building with a light-colored facade and many windows. The building has a sign that reads "PENTAX SURGICAL" on its side. A street lamp is visible in the foreground on the right.

About Us!

ABOUT US

WELCOME TO PENTAX SURGICAL

Established in 2008, PENTAX SURGICAL is on a mission to redefine the industry by proving that quality doesn't have to be sacrificed for competitiveness. Our commitment revolves around producing high-quality instruments at prices that keep you smiling. From our humble beginnings as a family-run veriure, we've evolved into a significant player in the market.

Our mission is clear to instill confidence in our clients by delivering Instruments of the highest quality at prices that make sense. We understand that the key to succes lies in customer satisfaction and loyalty, which forms the core of our vision. As we grow, our focus reamins on fostering strong relatioriships with our customers

At PENTAX SURGICAL, we firmly believe that quality is never accidental, it's the outcome of thoughtful and intelligent effort Each instrument we produce reflects this commitment to excellence, ensuring you receive products that exceed expectations. Join us on our journey as we continue to set new standards in quality and customer satisfaction.



PS-L-001

Handle Standard Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	20-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	20-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	20-200-40
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	20-200-41

Sets available in hard case or soft pouch



20-100-00

Mc - Intosh Blade Standard Illumination

Size: 0



Mc - Intosh Blade Standard Illumination

Size: 1



Mc - Intosh Blade Standard Illumination

Size: 2



Mc - Intosh Blade Standard Illumination

Size: 3



Mc - Intosh Blade Standard Illumination

Size: 4

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Mc – Intosh Laryngoscopes

Fiber Optic Illumination

Adult Size



Mc - Intosh Blade Fiber Optic Illumination
Size: 0



Mc - Intosh Blade Fiber Optic Illumination
Size: 1



Mc - Intosh Blade Fiber Optic Illumination
Size: 2



Mc - Intosh Blade Fiber Optic Illumination
Size: 3



Mc - Intosh Blade Fiber Optic Illumination
Size: 4

PS-L-002

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	21-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	21-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	21-200-40
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	21-200-41

Sets available in hard case or soft pouch

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



PS-L-003

Handle Standard Illumination

Size: Child / Infant

Handle options: Adult & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Child Handle	▪ 22-200-30
▪ Set of 3 Blades Size 0, 1, 2 and Child Handle	▪ 22-200-31
▪ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	▪ 22-200-40
▪ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	▪ 22-200-41

Sets available in hard case or soft pouch



Miller Blade Standard Illumination
Size: 00



Miller Blade Standard Illumination
Size: 0



Miller Blade Standard Illumination
Size: 1



Miller Blade Standard Illumination
Size: 2



Miller Blade Standard Illumination
Size: 3

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Miller Laryngoscopes

Fiber Optic Illumination

Adult Size



Miller Blade Fiber Optic Illumination
Size: 0



Miller Blade Fiber Optic Illumination
Size: 1



Miller Blade Fiber Optic Illumination
Size: 2



Miller Blade Fiber Optic Illumination
Size: 3



Miller Blade Fiber Optic Illumination
Size: 4

PS-L-004

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	23-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	23-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	23-200-40
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	23-200-41

Sets available in hard case or soft pouch

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Mc – Intosh Integrated Laryngoscopes

Fiber Optic Illumination

Adult Size



PS-L-005

Handle Fiber Optic Illumination

Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Adult Handle	▪ 24-200-30
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ 24-200-31
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ 24-200-40
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ 24-200-41

Sets available in hard case or soft pouch



Mc - Intosh Blade Fiber Optic Illumination
Size: 00



Mc - Intosh Blade Fiber Optic Illumination
Size: 0



Mc - Intosh Blade Fiber Optic Illumination
Size: 1



Mc - Intosh Blade Fiber Optic Illumination
Size: 2



Mc - Intosh Blade Fiber Optic Illumination
Size: 3



24-100-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

McCoy Flexi - Tip Laryngoscopes

Fiber Optic Illumination

Adult Size



PS-L-006

Handle Fiber Optic Illumination
Size: Adult

Standard Sets

Contents	Code
▪ Set of 2 Blades Size 2, 3 and Adult Handle	▪ 25-200-30
▪ Set of 2 Blades Size 3, 4 and Adult Handle	▪ 25-200-31
▪ Set of 3 Blades Size 2, 3, 4 and Adult Handle	▪ 25-200-40

Sets available in hard case or soft pouch



McCoy Blade Fiber Optic Illumination
Size: 1



McCoy Blade Fiber Optic Illumination
Size: 2



McCoy Blade Fiber Optic Illumination
Size: 3



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

McCoy Flexi - Tip Integrated Laryngoscopes

Fiber Optic Illumination

Adult Size



PS-L-007

Handle Fiber Optic Illumination
Size: Adult



McCoy Blade Fiber Optic Illumination
Size: 1



McCoy Blade Fiber Optic Illumination
Size: 2



McCoy Blade Fiber Optic Illumination
Size: 3



Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	26-200-30
Set of 2 Blades Size 3, 4 and Adult Handle	26-200-31
Set of 3 Blades Size 2, 3, 4 and Adult Handle	26-200-40

Sets available in hard case or soft pouch

Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



PS-L-008

Brass Otoscope Fiber Optic
Brass Otoscope Conventional



PS-L-009

ENT Set Complete



PS-L-010

Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available
- Fiber Optic Light Source



PS-L-011

Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



Color Options

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source

Standard Accessories

Nasal Speculum



PS-L-012



PS-L-013



PS-L-014

Bulbs



PS-L-015



PS-L-016



PS-L-017



PS-L-018

Ear Speculas



PS-L-019



PS-L-020



PS-L-021



PS-L-022

Ear Speculas & Testing Tube



PS-L-023



PS-L-024



PS-L-025



PS-L-026



PS-L-027



PS-L-028



PS-L-029



CONTACT US



+92 328 01 66 200



rafay@pentaxsurgical.com
pentaxsurgical@gmail.com



www.pentaxsurgical.com